

# Automatic Distillation Apparatus

## Product Brief

Smart D6 automatic intelligent distiller is a fully automatic finishing system designed according to the principle of automatic distillation and the national standard system of distillation. The system adopts a strong arm control core to ensure that the system has a strong control core to improve the overall corresponding speed, reliability and scalability; highly reliable weighing terminal detection technology and professional interference control technology to ensure the high reliability of terminal judgment; PC terminal and network platform can be accessed to improve the most convenient operation and test control technical title for you. The instrument has the functions of accurate temperature control, automatic distillation, circulating condensation, anti backdraft, explosion-proof boiling, automatic control of distillate volume, online cleaning of condensation pipeline, etc., and it can complete the distillation process automatically. It is suitable for distillation and digestion of volatile phenol, cyanide, ammonia nitrogen, COD and other items in the sample.



Serial number	project	Smart N6	Smart M6	Smart D6
One	temperature control	Room temperature to 0-270 °C (single hole adjustable)	Room temperature to 320 °C (single hole adjustable)	
Two	time control	0-120min (single hole adjustable)	0-999min (single hole adjustable)	
Three	Heating unit	Six	six	six
Four	power	Six joint total power 3000W (single hole 400W adjustable)		
Five	Heating up time	—	Within 15min	Within 15min
Six	Key control mode of distillation	—	Weighing terminal control	Weighing terminal control
Seven	Distillation end point accuracy	—	0.01 grams	0.01 grams
Eight	Distillation setting	—	1-600g adjustable	1-600g adjustable
Nine	Cooling water circulation	—	—	automatic
Ten	Receiving bottle specification	Capacity bottle of any specification		

## To configure

One set of automatic distiller, six 500L double neck round bottom flasks, six condensing tubes, six receiving bottles, one silica gel tube and one set of random tools.

## Performance characteristics

### Truly automatic

Automatic completion of heating, temperature control, cooling water refrigeration, cooling water cycle, distillation end point detection, without manual intervention

### Far infrared heating

The infrared ceramic mold is used to form heating components, which has the characteristics of high heat conversion efficiency, uniform heating, non explosive boiling, corrosion resistance, etc

### Accurate judgment of end point

The distillation end-point adopts the weighing end-point judgment technology, combined with the interference control technology independently developed, to ensure the reliability of the end-point judgment. The end-point can be adjusted arbitrarily, and the end-point accuracy can reach 0.01g

### Flexible control of distillation unit

The number of distillation units is six, single hole and single control, more than six groups of sample distillation can be completed at one time, and the heating temperature and time can be freely set

### Built in refrigeration unit

Built in compression refrigeration component, closed circulation of cooling water source, no need for external tap water, to meet the condensation work under various working conditions

### Dual mode adjustable

There are two working modes: weighing mode and time recording mode. One key switch is applicable to all kinds of distillation occasions. Each unit of working mode can be set independently, and different units can distill different samples

### One button peeling

Built in receiving bottle peeling function, one button peeling, one button start the experiment process, reduce the workload of the experimenter to no lower level

### Distillation liquid anti backdraft

The function of preventing backdraft of distillate can avoid the end of distillation and ensure the accuracy of test results.

## Applicable standard

HJ 503-2009	Water quality -- Determination of volatile phenol -- 4-aminoantipyrine spectrophotometric method
HJ 484-2009	Water quality -- Determination of chloride -- Volumetric and spectrophotometric methods
HJ 537-2009	Water quality -- Determination of ammonia nitrogen -- distillation neutralization titration method
HJ 745-2015	Determination of chloride and total chloride in soil spectrophotometric method



- It adopts 5.7-inch full touch 16 bit true color liquid crystal display and smart operation interface, which is easy to operate and truly realizes one key test