


# 8 Strip PCR Tubes Capping Machine

Axygen 8 strip pcr tubes



Strip tube channels



**Description**

The 8 Strip PCR Tubes Capping Machine design for pcr tube strip type close lid with strip.

**Parameter**

**Input power:** DC24/5A

**Mold:** Axygen 8 strip pcr tubes

**Channel:** 12

**Operate:**

1. Connect power supply
2. Turn on machine
3. Self detection, if error, will alarm itself
4. Put strips into channel
5. Place lids on top of strips
6. Press “**Start**” button 2 seconds to run

**Video illustration:**

<https://youtu.be/KwPOBV7djWk>