


Filters For Pipette Tip



High quality filters for 10/200/1000/5000/10000ul

No.	Specs	Package
1	For 10ul pipette tip (F10-16, F10-18)	1000Pcs/bag
2	For 200ul pipette tip (F200-45, F200-47)	1000Pcs/bag
3	For 1000ul pipette tip (F1000-7, F1000-72, F1000-74)	1000Pcs/bag
4	For 5ml pipette tip (F5000-90-40)	300Pcs/bag
5	For 10ml pipette tip (F10000-142-40)	300Pcs/bag

Cat. No.	Diameter (mm)	Thickness (mm)	Pore size (mm)	For tip size	Package
F10-14	1.4	3.5	10	10ul	10000Pcs/bag
F10-14s	1.4	1.6	10	10ul	10000Pcs/bag
F10-16	1.6	3.5	10	10ul	10000Pcs/bag
F10-18	1.8	3.5	10	10ul	10000Pcs/bag
F10-21	2.1	3.5	10	10ul	10000Pcs/bag
F10-22	2.2	3.5	10	10ul	10000Pcs/bag
F20-25	2.5	3.5	10	20ul	10000Pcs/bag
F10-25	2.5	3.5	10	10um	10000Pcs/bag
F10-26	2.6	3.5	10	10um	10000Pcs/bag
F50-34	3.4	3.5	10	50ul	10000Pcs/bag
F100-38	3.8	3.2	10	100ul	10000Pcs/bag
F100-38MY	3.8	3.2	10	10um	10000Pcs/bag
F43	4.3	3.5	10	10ul	10000Pcs/bag
F200-45	4.5	3.5	10	200ul	10000Pcs/bag
F200-47	4.7	3.5	10	200ul	10000Pcs/bag
F200-47MY	4.7	3.2	10	200ul	10000Pcs/bag
F200-49	4.9	3.5	10	200ul	10000Pcs/bag
F200-51	5.1	3.5	10	10um	10000Pcs/bag
F55	5.5	3.5	10	10um	10000Pcs/bag
F58	5.8	3.5	10	10um	10000Pcs/bag
F500-60	6.0	4.0	10	500ul	10000Pcs/bag
F1000-70	7.0	4.0	10	1000ul	10000Pcs/bag
F1000-71	7.1	4.0	10	1000ul	10000Pcs/bag
F1000-72	7.2	4.0	10	1000ul	10000Pcs/bag
F1000-74	7.4	4.0	10	1000ul	10000Pcs/bag
F78-80	7.8	8.0	10	10ul	10000Pcs/bag