


Immunoassay Analyzer



Overview



Features:

The lowest detection limit: 0.05% of the full-scale signal value

Instrument stability: 10 measurements, TC ratio CV less than 0.5%

Instrument SNR: 0.01% (1:10000)

Fluorescence detection bands: 365/610, 470/540, 525/610, 610/690 (nm)

Colloidal gold light source band: 525 (nm)

Support network connection: WIFI, GPRS

- WIFI connection: LCD screen select WIFI hotspot, enter password
- GPRS connection: configure the flow card, automatically connect to the 4G base station
- Automatic upload of detection data
- Automatic fault alarm, remote diagnosis
- Automatic download of standard curve

Parameters

Dimension	220*100*75mm
Weight	< 1kg
Touch screen	4.3 in
Resolution	480*272 px
Power	6000mAh lithium
Working hour	continuously working 10 hours
Connection port	USD, WIFI, GPRS, BlueTooth
Running environment	5°C~40°C , humidity 10~80%
Memory	8,000 records
Working model	Single machine, wifi connect with PC software
Program	1200 programs & lots
Cassette	1 cassette
T lines	Maximum 4 T lines
Calibration	Internal calibration

Components

Item	Description	Number
Housing	Portable, plastic	1
Detection model	Fluorescence/colloidal gold	1
Touch screen	4.3 in	1
Memory	Storage records	1
PCB	Control board	1
Communication module	wifi	1
Communication module	GPRS	1
Communication module	Bluetooth	1
Printer	External wifi printer	1
Printer	External serial port printer	1
Barcode	HM code scanner	1
Data cable	USB cable	1
Battery	Lithium cell 18650, 6000mAh	2
Power cable	12V, 2A DC adapter	1