

Nucleic Acid Purification Spin Columns

Antiteck nucleic acid purification spin columns are based on Antiteck's proprietary spin column-based silica membrane purification technology. In addition to the spin column itself, each spin column comes with one collection tube, silica membrane and fix ring.



Ordering Information:

Cat.#	Description	Volume	Elution volume	Qty.
RP20-A-N	RNA Extraction Columns, Capped spin columns, Green fixing rings	2 mL, 800 μ L	\sim 20 μ g	500 Pcs/PK
NP20-A	Genomic DNA Extraction Columns, Capped spin columns, Green fixing rings	2 mL, 800 μ L	\sim 20 μ g	500 Pcs/PK
NP20	Genomic DNA Extraction Columns, Capless spin columns, Green fixing rings	2 mL, 800 μ L	\sim 20 μ g	500 Pcs/PK
NP30	Plasmid Miniprep Columns, Capless spin columns, Blue fixing rings	2 mL, 800 μ L	\sim 30 μ g	500 Pcs/PK
NP30-A	Plasmid Miniprep Columns, Capped spin columns, Blue fixing rings	2 mL, 800 μ L	\sim 30 μ g	500 Pcs/PK
HP50	Plasmid Miniprep Columns, Large-scale, Capless spin columns, White fixing rings	2 mL, 800 μ L	\sim 80 μ g	500 Pcs/PK
HP50-A	Plasmid Miniprep Columns, Large-scale, Capped spin columns, White fixing rings	2 mL, 800 μ L	\sim 80 μ g	500 Pcs/PK
NP05	DNA Cleanup & Gel Purification Columns, Mini spin columns, Capless spin columns, White fixing rings	2 mL, 800 μ L	\sim 10 μ g	500 Pcs/PK
NP10	DNA Cleanup & Gel Purification Columns, Capless spin columns, Yellow fixing rings	2 mL, 800 μ L	\sim 10 μ g	500 Pcs/PK
NP10-A	DNA Cleanup & Gel Purification Columns, Capped spin columns, Yellow fixing rings	2 mL, 800 μ L	\sim 10 μ g	500 Pcs/PK
MP05	Ultra micro DNA extraction columns, Capped spin columns, White fixing ring.	2ml, 800 μ L	\sim 5 μ g	500 Pcs/PK
NP100	Plasmid Midiprep Columns	15 mL, 4 mL	\sim 100 μ g	50 Pcs/PK
HP200	High-yield Plasmid Midiprep Columns	15 mL, 4 mL	\sim 200 μ g	50 Pcs/PK
NP500	Plasmid Maxiprep Columns	50 mL, 18 mL	\sim 500 μ g	20 Pcs/PK
HP1000-P	High-yield Plasmid Maxiprep Columns	50 mL, 18 mL	\sim 1 mg	20 Pcs/PK
RS30-2	Plasmid Megaprep Columns	30 mL	500-1000 μ g	20 pcs/Box
RS30-22	Plasmid Gigaprep Columns	300 mL	5-10 mg	10 Pcs/PK
EP50-06	Large-Volume Plasmid Miniprep Columns, 6 mL extension tubes, capless spin columns, white fixing rings	\sim 80 μ g	50/Box	50 Pcs/PK
EP50-06A	Large-Volume Plasmid Miniprep Columns, 6 mL extension tubes, capped spin columns, white fixing rings	\sim 80 μ g	50/Box	50 Pcs/PK
EP50-12	Large-Volume Plasmid Miniprep Columns, 12 mL extension tubes, capless spin columns, white fixing rings	\sim 120 μ g	50/Box	50 Pcs/PK
EP50-12A	Large-Volume Plasmid Miniprep Columns, 12 mL extension tubes, capped spin columns, white fixing rings	\sim 120 μ g	50/Box	50 Pcs/PK
DNAK9602-N	96-Well Extraction Plates, Semi-Skirted	1.0 mL/well	\sim 15 μ g/well	4 Pcs/PK
DNAK9607-N	96-well Extraction Plates, Full-Skirted	1.5 mL/well	\sim 15 μ g/well	4 Pcs/PK

ANTITECK LIFE SCIENCES LIMITED

Add: A1-519, XingGang GuoJi, Yingbin Road, Huadu, Guangzhou, China

Tel: +86 133 1286 3972

Web: <https://antiteck.com>

E-mail: info@antiteck.com